FastFinder

Automated PCR analysis

Instructions for use

Software release version 3.6 (released June 2020)
IFU version 1.0

Kempische Steenweg 303/105, 3500 Hasselt, Belgium
Office: +32 78 48 44 79
Support: +32 78 48 46 55
support@ugentec.com

* Due to regulatory reasons, this software cannot be used for diagnostic purposes outside of the European Union.

1 Intended Purpose

FastFinder Platform, as part of the FastFinder System, is intended to be used by trained laboratory professionals as decision-support software for the analysis and interpretation of PCR instrument results. The FastFinder Platform hosts, processes and displays the results of assay-specific applications (assay plugins), which contain assay specific algorithms and decision mechanisms that convert PCR raw data into final test results with minimal manual intervention.

It is assumed that the user is a professional with the necessary clinical and/or laboratory expertise, background and experience in molecular biology that has received adequate training on the FastFinder platform.

All final results should be interpreted in conjunction with the patient's clinical signs and symptoms, medical history and other laboratory data available to the physician.

2 FastFinder Platform; an introduction

FastFinder is a platform on which UgenTec hosts assay-specific “plugin” applications that support the analysis and interpretation of PCR instrument results. This platform does not generate end-results as stand-alone software. Only in combination with a PCR Assay and an Assay Plugin, this generic FastFinder Platform can turn into a system for automating PCR analysis interpretation for a specific assay.

An Assay Plugin is a self-contained file, that contains an assay-specific algorithm and decision mechanisms which allow to convert the raw PCR data into end results. Each Assay Plugin is designed and developed in collaboration with the assay developer. The algorithm is the data analysis component of the Assay Plugins, which contains a model, that is trained via machine learning technology, to differentiate positive curves from negative curves and which determines the C_q values of the PCR curves. The decision tree is the rules engine component of the Assay Plugin which determines how the Algorithm results are to be interpreted by the FastFinder System. It contains criteria to determine a final result as defined in the assay instructions for use or standard operating procedures. So, the result is an analysis workflow that begins with the PCR raw data file and ends with a sample result.


The FastFinder Platform hosts Assay Plugins for clinical *in vitro* diagnostic use as well as for research use only (RUO) or Assay Plugins for molecular testing in the veterinary, agricultural or food industry. The Product Specification Sheet of the Assay Plugin clearly contains a statement to indicate the intended use of that particular Assay Plugin.

The FastFinder Platform is designed to feel as familiar and intuitive as an internet browser with all the samples and analyses displayed as hyperlinks accessible with a single click. From the homepage, all available or ongoing analyses can be accessed and in the uppermost pane there are separate FastFinder Platform modules devoted to *Analyses*,

Devices, Assays, QC, Archive, Settings and *Help*. Multiple tabs enabling multiple runs to be analysed and compared can be opened at the same time.

FastFinder has a strong security model implementing role-based user access control and data access segregation, with supporting audit trails.

The following symbols are used in FastFinder:

	Edit function
	Select
	Zoom
	Approved
	Authorised or Default
	Restricted
	Drop down menu
	Control assay
	Research Use Only assay
	Resolve item(s)
	Audit trail
	Reject sample
	Approve sample
	To be reviewed
	To be repeated
	Date of manufacture

	Manufacturer
	Serial number
	Reset
	Set filters
	Troubleshooting
	Search
	One click analysis

Key concepts:

- **Assay Plugin:** Self-contained assay-specific application that is part of the FastFinder System and contains assay-specific algorithms and decision mechanisms which allow to convert the raw PCR data into end results.
- **Algorithm:** The data analysis component of an Assay Plugin, which contains a model, that is trained via machine learning technology, to differentiate positive curves from negative curves and which determines the Cq values of the PCR curves.
- **Decision tree:** The decision tree is the rules engine component of the Assay Plugin which determines how the Algorithm results are to be interpreted by the FastFinder System. It contains criteria to determine a final result as defined in the assay instructions for use.
- **Laboratory environment:** the shared workspace for all FastFinder users in your laboratory, pertaining to a single database and a set of users managed centrally from the admin module (see section 2.2 Managing users and permissions in the technical note).

3 Minimum system requirements

Hardware requirements

- Processor: 2 GHz, 2 GB RAM
- Disk space: 10 GB
- Internet Connection Cable or DSL
- Screen resolution: min. 1366 x 768 pixels, max. 1920 x 1080 pixels

Operating system

- Windows 7 or higher

Supported browsers

The FastFinder Admin, where users are managed, requires:

- Internet Explorer 11 or newer
- Microsoft Edge 25 or newer
- Firefox 45 or newer
- Google Chrome 47 or newer


FastFinder can be installed by downloading the installer in the registration email and open the executable file and following the installation wizard. When the user is able to login, perform an analysis and create a report, the installation was successful.

4 Supported PCR cyclers and files

Device	Filetype(s)
Applied Biosystems ABI® 7500 (fast) (Dx)	.sds (software v1.x) .eds (software v2.x)
Applied Biosystems ViiA7	.eds
Applied Biosystems QuantStudio 3 (seen as QS5) and QuantStudio 5	.eds
Applied Biosystems QuantStudio 5 Dx	.eds
Applied Biosystems QuantStudio 6 and 7	.eds
Applied Biosystems QuantStudio 7 Pro	.eds
Applied Biosystems QuantStudio Dx	.eds
Abbott® m2000	.sds
Agilent® Mx3000P	.mxd
Agilent® Mx3005P	.mxd
Cobas® z 480	.ixl, .xml
BioRad® CFX (CFX 96, CFX 384)	.pcrd
Roche LightCycler® 480 Type I & II	.ixl, .xml
Qiagen Rotor-gene 6000 & Q	.rex
Cepheid SmartCycler®	.csv*
BMS micPCR	.micrun
FastTrack Cycler	.ftcrun
RIDA®CYCLER	.rcyclerrun
AriaDx	.amxd
AriaMx	.amxd
LightCycler 96	.lc96p

QuantStudio 12K Flex	.eds
IntelliQube	.csv

* Ensure that SmartCycler CSV files are exported with at least the following two options checked:

5 Incidents and Support

For technical support in setting up and operating your FastFinder laboratory environment or in case of malfunctioning, please visit our support center <https://support.ugentec.com/> or contact support@ugentec.com. For a more detailed explanation of all FastFinder Platform functionalities, please consult the technical note (see help module of the FastFinder Platform).

Any incidents with the FastFinder Platform need to be reported to our support center <https://support.ugentec.com/> or support@ugentec.com. Serious incidents in relation to the FastFinder Platform need to be reported to the competent authority of the country in which the end user and/or patient is established.

6 Product classification and certification

FastFinder Platform is created under a certified ISO 13485 Quality Management System.

EEA (European Economic Area)

FastFinder Platform is marketed in Europe as a general in vitro diagnostic medical device in conformity with the essential requirements and provisions of the Council Directive 98/79/EC concerning in vitro diagnostic medical devices.

USA (United States of America)

Based on definitions set forth in the FDA "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices", the FastFinder Platform is marketed in the USA as an in vitro diagnostic medical device system component, with a minor level of concern and subject to FDA general controls and software quality system regulations for medical devices. The level of concern of the Assay Plugins that are hosted on the FastFinder Platform are determined by the device classification of the corresponding in vitro diagnostic assay.

Other regions and countries

In any country or region not covered above, FastFinder Platform is marketed for Research Use Only.

7 Product functionalities

7.1 Analysis types

FastFinder Platform supports the following PCR analysis types:

- amplification curve analysis
- melting curve analysis
- absolute quantification analysis
- relative quantification analysis

All assay specific information (device type, assay type, used channels and targets, algorithm and decision tree) is fixed of the Assay Plugin. The FastFinder Platform only visualises the results of the algorithm and decision tree as defined in the Assay Plugin. FastFinder calculates concentrations and relative quantifications based on the information in the Assay Plugin.

For some device types and/or analysis types, additional calibration information is necessary in order to get the correct interpretation of the results. The customer is responsible to link the correct calibration files and for the quality of the calibration experiments. Also for proper lot tracking, the correct lot needs to be selected for each assay.

FastFinder allows default settings for color compensation, T_m calibration and used lots. The user is responsible for the correct use of this functionality.

Errors or warnings encounters by the FastFinder Platform are displayed in the user interface. Warnings are indicated by an orange symbol, errors by a red symbol. It is the sole responsibility of the user to proceed with an analysis for which warnings raised. In case of errors, the user cannot proceed in FastFinder.

7.2 Assay Plugins

An Assay Plugin is the digital representation of a physical laboratory test or commercial kit. It contains all assay specific information that is needed for the analysis and interpretation of a PCR assay:

- creator
- device type
- analysis type
- used channels and targets
- sample types and controls
- replicates
- algorithm
- decision tree
- custom parameters and calculations
- FastFinder restrictions

An assay plugin is a self-contained piece of software. They are encrypted. This means that they can never be changed once they have been created. The FastFinder Platform can host IVD as well as RUO Assay Plugins. Only in combination with an IVD Assay and Assay Plugin, this generic FastFinder Platform can turn into an IVD system for automating PCR analysis and interpretation for a specific IVD assay. Automatic assignment of Assay Plugins is possible by the use of nametags. Incorrect settings or use of nametags can result in a wrong plate assignment. The user is responsible for the correct use of nametags

Backwards compatibility

UgenTec can confirm that Assay Plugins, released on previous FastFinder Platform versions, will produce the same end result on this version of the FastFinder platform. A FastFinder Platform update does not affect the safety, effectiveness, (clinical) functionality, or performance specifications that are directly associated with the intended use of the Assay Plugins that are hosted on FastFinder Platform.

7.3 QC module

The optional QC module tracks selected control sample types/targets in time. The FastFinder Platform makes use of the Westgard rules in combination with a fix or floating mean/SD calculation. For more information about the Westgard rules, please visit <https://www.westgard.com/>. Configuration for the QC follow up can be changed during lot creation. The QC tracking is scoped per device/device type - assay - lot - target combination.

The FastFinder Platform only considers QC results that were analysed in FastFinder. For a correct use of the QC module, the user is responsible to analyse the PCR samples in a chronological order.

7.4 2-step validation

FastFinder platform supports two-step validation (optional). In the settings, the user rights to be a first or second reviewer can be assigned to a specific function within your laboratory environment. The second reviewer can review all results, but gets an additional overview from all the changes and remarks made by the first reviewer.

7.5 Exports

FastFinder Platform supports a standard pdf, csv and ASTM export. Custom implementations can be created upon request (please contact support@ugentec.com for more information about the different options). Pdf exports can be generated at any time in an analysis, csv exports are only allowed for authorised analyses.

The standard pdf and csv exports are verified and validated as part of the FastFinder platform. All custom implementations need to be validated by the end user before clinical use. The ASTM export outputs a defined set of information, and the user's LIMS system need to filter the relevant data for your application. It is the end user's responsibility to validate the LIMS connection before clinical use. Technical instructions for the ASTM export are available on request via support@ugentec.com.

Exports can be modified manually after generation outside of FastFinder and hence do not necessarily represent the actual state of the data in the FastFinder Platform database at time of generation. The content of exports is the sole responsibility of the user.

7.6 Archive

All analyses are stored in the archive module. Even rejected analyses can be consulted here. Authorised analyses are read-only and cannot be edited anymore. Analyses that are still in progress cannot be edited by two users at the same time.

7.7 Custom settings

All settings in FastFinder are on a laboratory level, so applicable for all users in your lab. Following settings can be customised:

- User authorisations
- General lab settings
- Default graph and results overview settings
- File watcher
- Roche FLOW integration
- QC settings (applicable Westgard rules, optional)
- LIMS settings (optional)
- Network configuration (proxy settings)
- Report settings

For more information about these settings, please consult the technical note.

If the IFU of the assay that is analysed with the FastFinder system contains instructions on custom settings that are needed in order to correctly analyse and interpret that particular assay, it is the user's responsibility to comply with the instructions provided in the assay's instructions for use.

8 Performance Characteristics

The FastFinder Platform cannot demonstrate clinical effectiveness as standalone software. Analytical and clinical performance is determined by the associated IVD assays of the Assay Plugins, which are hosted on the FastFinder Platform. Analytical and clinical performance of these IVD assays is described in their respective instructions for use.

9 Document versions

Version	Description	Date
1.0	New document	June 1, 2020

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